


StyleView Drawer Reference Guide





ergotron


Standard Configurations





Drawer Sizes




























Note:

Reconfigurable divider inserts are included with all drawers; inserts are interchangeable









Total weight capacity per drawer. ≤0.25 lbs (0.1 kg)



Configure Your Own

- Ergotron offers a range of drawer accessories, which can be purchased separately and added to any standard cart configuration
- All StyleView drawer systems consist of one "primary" drawer followed by additional "supplemental" drawers
- Choose between auto-locking and manual key-locking drawer options



Configure Your Own

- Any combination of drawers is allowed, provided they do not exceed the equivalent of 5 single-height rows
- Double-height drawers (1x1 Tall and 2x1 Tall) equal two single-height rows
- 3 or more single-height rows will reduce the cart's height-adjustment range
- A Travel-Stop accessory is recommended when 3 or more rows are installed, preventing drawer contact with the cart's base



Get More Info:

If your application requires a non-standard drawer layout, contact an Ergotron representative, who can help you configure your own system with compatible accessory options.

